

PneumoDart

CE REF TM-317

NSN: 6515-01-655-9514

Vol. 800 per case, Total Weight: 16kg/35.2lbs


The PneumoDart is a reliable solution for treating a tension pneumothorax injury, or sucking chest wound. Our robust innovation helps you to get it right, first time.


It offers better efficiency, safety and important physical features. Typical methods for treating this kind of injury include the application of a plastic seal over the chest wound to inserting a needle catheter into the chest to drain air out of the chest and relieve pressure. However, successful application of a traditional pneumothorax needle has proven difficult.

Our needle is 3.25" which is the optimal length as recommended by the Tactical Combat Casualty Care (TCCC) who inspired the PneumoDart and we developed it with the US Department of Defense during a robust testing and evaluation phase.

Details:

- 14ga x 3.25"
- Veress safety tip to reduce the risk of puncturing internal organs when inserted
- Visible cap head to monitor air or blood leaving the chest
- Luer lock cap allows for an active drawing or elimination of air or blood from the chest
- Seven year shelf life





Product Range Overview

Product Name	Part No.	NSN	CE Mark	Volume	Weight
4"Trauma Bandage	TM-204	6510-01-558-4108	Yes	100 per case	8kg/17.6lbs
6"Trauma Bandage	TM-206	6510-01-558-4114	Yes	50 per case	6kg/13.2lbs
ABD Trauma Bandage	TM-212	6510-01-541-8121	Yes	40 per case	7kg/15.4lbs
EZ Gauze	TM-215	6510-01-605-2234	Yes	200 per case	9kg/19.8lbs
TPAK	TM-303	6515-01-541-0635	Yes	1000 per case	15.5kg/34lbs
PneumoDart	TM-317	6515-01-655-9514	Yes	800 per case	16kg/35.2lbs
Xtricate	TM-511	6530-01-477-8525	Yes	1 per case	8kg/17.6lbs
TVAC	TM-909	6515-01-540-7206	Yes	20 per case	4.8kg/10.5lbs
Piranha Trauma Shears Large (7.25")	TM-707	6515-01-538-9276	Yes	20 per case	1.6kg/3.5lbs
Piranha Trauma Shears Small (5.50")	TM-705	6515-01-538-9276	Yes	30 per case	1.4kg/3lbs